



ACCESSORIES POLICY

The Budapest Home Office can only support the safety and effectiveness of the limb harness (LH) and head harness (HH) that is part of the registered device as made by us.

The **only** accessories that are registered as part of the original device you can see below.

SCIO HEAD AND LIMB HARNESS

The head and limb harnesses are shipped directly to the customer from the Budapest Home Office. These are marked with the batch number LHxxxxxx and HHxxxxxx. They look like the head and limb straps below:





The use of any other accessory including head or limb harness is at the risk of the user. Please be aware that only the above pictured head and limb harnesses are tested for safety and effectiveness and are registered with the SCIO.

Please note that the Budapest Home Office cannot take any responsibilities for use of any other harnesses and/or accessories made by a Third Party. Using such harnesses and/or accessories are at the user's own risk.

Please always note that the Third Party is responsible for getting their harnesses and/or accessories safety tested by a proper testing body.

The user has the right to ask the Third Party for a safety certificate of the harnesses and/or accessories. Please note that as of 1 November 2008, this policy comes into effect. However, Third Parties are given six (6) months to get the safety certificate. Therefore, from 1 May 2009 the user should be asking for a safety certificate of harnesses and/or accessories.

NEW SCIO LIMB HARNESS

The Budapest Home Office is proud to present the new and more improved limb harness for your device, designed to be more comfortable and flexible to use, while maintaining a high-quality level. This harness is unique on the market today, not only because it offers you simpler ways of working with a patient, but because it uses the highest quality in materials and further maintenance.

In our 20 years of experience in the field of energetic medicine and biofeedback, volcanic rubber has proved to be the safest and the most conductive material for the harness, therefore, we have kept the same material in manufacturing the new harnesses.






With the covering of the metal part that at the previous model was in direct contact with the skin, we have created an even safer harness, eliminating the risk of bruises, sore or blisters when working with patients with electrical hyper reactivity.


The new design of the harness will make it more comfortable for you to deal with the patient during the biofeedback session. The cables that connect the plug to the straps are now made from a flexible and more resistant material. Also, you will notice that the bananas are painted the four colors that you have become very familiar with, red, blue, yellow and black, colors that indicate which extremity to connect the strap to. In the past, that was indicated by the small colored buckle located on the straps. In order to prevent tangling of the harness and provide easy handling, the cables have been tied together at an appropriate length, so you can connect your patient with no effort. The length can be adjusted according to your needs.

The connection of the straps to the cables is now made through banana plugs. This way, if during the session your patient needs to be disconnected from the harness for any reason, you can do so with a simple click.






You can either remove the strap or unplug it, as seen in the pictures below.



When your patient returns and is ready to go back on the device, you simply re-connect the patient in the correct way , following the colored hints.

Another exciting advantage of this design is the fact that it increases the duration of the harness' life cycle. If, for any reason, your harness is damaged, you no longer need to exchange the whole harness and purchase a new one. The new design allows you to order only the damaged part and replace it, the cost of maintenance being significantly reduced.

We used the same serial plug for the new harness, so it is compatible with the SCIO device.



We hope that with the new model of the limb harness, we have improved the quality of your work with the SCIO Universal Electrophysiological Biofeedback System, so that you and your patients will be able to enjoy these advantages and have a more pleasant biofeedback experience. We are constantly striving to offer you the highest quality in services and products and bring you tomorrow's medicine...today!



EDUCTOR ACCESSORIES


The **only** accessories that are registered as part of the original device is the head and limb harnesses and the test plate, which are shipped directly to the customer from the Budapest Home Office. These are marked with the batch number LHxxxxxx and HHxxxxxx and TPxxxxxx. They are made from the same material as the harnesses of the SCIO and look like the head and limb straps below:



EDUCTOR Head Harness



EDUCTOR Limb Harness



EDUCTOR Test Plate